

DGK-α (A340) polyclonal antibody

Catalog: BS1960

Host: Rabbit

Reactivity: Human, Rat

BackGround:

Diacylglycerol kinases (DGKs) phosphorylate diacylglycerol (DAG) to produce phosphatidic acid. DAG and phosphatidic acid are lipids that act as second messengers in signaling cascades. DGK-alpha influences cell activation and secretion of lethal exosomes, which in turn control cell death. DGK-beta is abundant in restricted brain regions such as the caudate putamen and olfactory tubercle. DGK-gamma encodes full-length and truncated transcripts that are present in a range of human tissues, with greatest expression observed in retina. DGK-delta is most abundant in skeletal muscle. DGK-epsilon shows specificity for arachidonyl-containing diacylglycerol and is expressed predominantly in testis. DGK-zeta is most abundant in brain and muscle. DGK-eta is closely related to DGK-delta. DGK-theta is most abundant in the cerebellum and hippocampus.

Product:

Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2

Molecular Weight:

~ 90 kDa

Swiss-Prot:

P23743

Purification&Purity:

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).

Applications:

WB: 1:500~1:1000

IHC: 1:50~1:200

Storage&Stability:


Store at 4 °C short term. Aliquot and store at -20 °C long

term. Avoid freeze-thaw cycles.


Specificity:

DGK- α (A340) polyclonal antibody detects endogenous levels of DGK- α protein.

DATA:



Western blot (WB) analysis of DGK-α (A340) pAb at 1:2000 dilution Lane1:C6 whole cell lysate(40ug) Lane2:U-87MG whole cell lysate(40ug) Lane3:Myla2059 whole cell lysate(40ug) Lane4:HuT78 whole cell lysate(10ug)



Immunohistochemistry (IHC) analyzes of DGK- α (A340) pAb in paraffin-embedded human tonsil carcinoma tissue at 1:50.

Note:

For research use only, not for use in diagnostic procedure.

Bioworld Technology, Inc.

 Add:
 1660 South Highway 100, Suite 500 St. Louis Park, MN 55416,USA.

 Email:
 info@bioworlde.com

 Tel:
 6123263284

 Fax:
 6122933841

Bioworld technology, co. Ltd. Add: No 9, weidi road Qixia District Nanjing, 210046,

P. R. China. Email: <u>info@biogot.com</u> Tel: 0086-025-68037686 Fax: 0086-025-68035151