# PRODUCT DATA SHEET



Bioworld Technology,Inc.

# Recombinant G-CSF, Rat (HEK 293-expressed)

Catalog Number: BK0305-10µg Source: HEK 293 Quantity: 10µg

### **Description:**

Among the family of colony-stimulating factors, Granulocyte Colony-Stimulating Factor (G-CSF) is the most potent inducer of terminal differentiation of leukemic myeloid cell lines into granulocytes and macrophages. G-CSF synthesis can be induced by bacterial endotoxins, TNF, Interleukin-1 and GM-CSF. Prostaglandin E2 inhibits G-CSF synthesis. In epithelial, endothelial, and fibroblastic cells, secretion of G-CSF is induced by Interleukin-17.

#### **Molecular Weight:**

25~28 kDa, observed by reducing SDS-PAGE.

#### **Purity:**

> 95% as analyzed by SDS-PAGE and HPLC.

# **Biological Activity:**

ED50 < 5 pg/ml, measured in a cell proliferation assay using NFS-60 cells.

# **Physical Appearance:**

Sterile Filtered White lyophilized (freeze-dried) powder.

#### **Formulation:**

Lyophilized after extensive dialysis against PBS.

#### **AA Sequence:**

IPLLTVSSLPPSLPLPRSFLLKSLEQVRKIQARN-TELLEQLCATYKLCHPEELVLFGHSLGIP-KASLSSCSSQAL-QQTKCLSQLHSGLFLYQGLLQALAGIS-SELAPTLDMLHLDVDNFATTIWQQMES-LGVAPTVQPTQSTMPIFTSAFQRRAGGVLV-TSYLQSFLETAHHALHHLPRPAQKHFPESLFISI

#### **Endotoxin:**

< 0.2 EU/µg, determined by LAL method.

#### **Reconstitution:**

Reconstituted in ddH2O or PBS at 100 µg/ml.

## **Storage:**

Lyophilized recombinant Rat Granulocyte Colony-Stimulating Factor (G-CSF) remains stable up to 6 months at -80 °C from date of receipt. Upon reconstitution, Rat Granulocyte Colony-Stimulating Factor (G-CSF) should be stable up to 1 week at 4 °C or up to 2 months at -20 °C.

# **Usage:**

This material is offered by USA Bioworld biotech for research, laboratory or further evaluation purposes. For research use only.