

## PRODUCT DATA SHEET

Bioworld Technology CO., Ltd.



### TIMP-4 (S206) Peptide

Cat No.: BS1477P

#### Background

TIMP-1, TIMP-2, TIMP-3 and TIMP-4 (for tissue inhibitor of metalloproteinases-1, -2, -3 and -4) complex with metalloproteinases such as collagenases, gelatinases and stromelysins, resulting in irreversible inactivation of the metalloproteinase. TIMP-1 was found to be identical to EPA (erythroidpotentiation activity). Parathyroid hormone has been shown to be a regulator of TIMP-2 in osteoblastic cells. TIMP-3 may be involved in regulating trophoblastic invasion of the uterus as well as in regulating remodeling of the extracellular matrix during the folding of epithelia, and in the formation, branching and expansion of epithelial tubes. TIMP-4 is most highly expressed in heart and low levels of TIMP-4 are expressed in liver, brain, lung, thymus and spleen.

#### Swiss-Prot

Q99727

#### Applications

Blocking

#### Specificity

This peptide can be used with studies using BS1477 TIMP-4 (S206) pAb.

#### Purification & Purity

Synthetic peptide TIMP-4 (S206). (Note: the amino acid sequence is proprietary). The purity is > 98%.

#### Product

1 mg/ml in DI water.

#### Storage & Stability

Store at 4 °C short term. Aliquot and store at -20 °C long term. Avoid freeze-thaw cycles.

#### Research Use

For research use only, not for use in diagnostic procedure.

Bioworld Technology, Inc.  
1660 South Highway 100, Suite 500 St. Louis Park, MN  
55416, USA. Email: [info@bioworld.com](mailto:info@bioworld.com)  
Tel: 6123263284 Fax: 6122933841

Bioworld technology, co, Ltd.  
No 9, weidi road Qixia District Nanjing, 210046,  
P, R.China. Email: [info@biogot.com](mailto:info@biogot.com)  
Tel: +86-025-68037686 Fax: +86-025-68035151