

**ENaC δ (P446) polyclonal antibody**

Catalog: BS3070

Host: Rabbit

Reactivity: Human

BackGround:

The epithelial sodium channel (ENaC) is a member of the ENaC/DEG superfamily that is located on the apical surface of cells. ENaC mediates sodium reabsorption in kidney, distal colon, lung, ducts of exocrine glands and other organs. ENaC is formed by heteromultimerization of four homologous subunits, α , β , γ and δ . The most frequently formed heterotetramer consists of 2 α , 1 β , and 1 γ subunit, but the α subunit can be replaced by a δ subunit. The α ENaC gene maps to human chromosome 12p13, and expresses a glycosylated protein. Both the β and γ ENaC genes map to human chromosome 16p12, and the γ ENaC transcript is detected as a glycosylated protein. The carboxy-terminus of all ENaC subunits contains PY motifs, which interact with the ubiquitin protein ligase, Nedd4, to regulate intracellular sodium concentrations. Gain-of-function mutations involving the PY motif cause Liddle's syndrome, an autosomal dominant form of hypertension, resulting from excessive renal sodium absorption. Conversely, ENaC loss-of-function mutations result in pseudohypoaldosteronism type I, a disorder characterized by salt wasting and hypotension.

Product:

Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2

Molecular Weight:

~80 kDa

Swiss-Prot:

P51172

Purification&Purity:

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).

Applications:

IHC: 1:50~1:200

IF: 1:50~1:200

Storage&Stability:

Store at 4 °C short term. Aliquot and store at -20 °C long term. Avoid freeze-thaw cycles.

Specificity:

ENaC δ (P446) polyclonal antibody detects endogenous levels of ENaC δ protein.

DATA:**Note:**

For research use only, not for use in diagnostic procedure.

Bioworld Technology, Inc.

Add: 1660 South Highway 100, Suite 500 St. Louis Park, MN 55416, USA.

Email: info@bioworld.com

Tel: 6123263284

Fax: 6122933841

Bioworld technology, co. Ltd.

Add: No 9, weidi road Qixia District Nanjing, 210046, P. R. China.

Email: info@biogot.com

Tel: 0086-025-68037686

Fax: 0086-025-68035151