

## PRODUCT DATA SHEET

Bioworld Technology,Inc.

# **Ghrelin polyclonal antibody**

Catalog: BS6800 Host: Rabbit Reactivity: Human, Rat

#### **BackGround:**

Small synthetic molecules called growth-hormone secretagogues (GHSs) act through GSH-R to stimulate the release of GH from the pituitary. Ghrelin is an endogenous levels of ligand for the growth hormone secretagogue receptor (GHS-R). Ghrelin and GHRH are involved in the regulation of GH release from the pituitary. GHRH exerts its action through high-affinity GHRH receptors (GHRH-R) present in the anterior pituitary. The acylated peptide of ghrelin specifically releases GH both in vivo and in vitro and is found in stomach tissue. GH plays a crucial role in stimulating and controlling the growth, metabolism and differentiation of many mammalian cell types by modulating the synthesis of multiple mRNA species. The gene encoding ghrelin maps to human chromosome 3p26-p25.

#### **Product:**

Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2

## **Molecular Weight:**

~ 13 kDa

## **Swiss-Prot:**

O9UBU3

## **Purification&Purity:**

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).

#### **Applications:**

WB: 1:500~1:2000 IHC: 1:50~1:200


## **Storage&Stability:**

Store at  $4 \,\mathrm{C}$  short term. Aliquot and store at  $-20 \,\mathrm{C}$  long term. Avoid freeze-thaw cycles.

#### **Specificity:**

Ghrelin polyclonal antibody detects endogenous levels of Ghrelin protein.

#### **DATA:**



Western blot analysis of extracts of HeLa cell lines, using Ghrelin antibody.

## Note:

For research use only, not for use in diagnostic procedure.

Bioworld Technology, Inc.

Add: 1660 South Highway 100, Suite 500 St. Louis Park,

MN 55416,USA.

Email: <u>info@bioworlde.com</u>

Tel: 6123263284 Fax: 6122933841 Bioworld technology, co. Ltd.

Add: No 9, weidi road Qixia District Nanjing, 210046,

P. R. China.

Email: <u>info@biogot.com</u>
Tel: 0086-025-68037686
Fax: 0086-025-68035151