

# NUP153 polyclonal antibody

Catalog: BS90977

Host: R

Rabbit

Reactivity: Human

# **BackGround:**

Nuclear pore complexes (NPCs) are the channels for the bi-directional movement of macromolecules between the nucleus and cytoplasm, and contain more than 100 different subunits. Many of them belong to a family called nucleoporins, which are characterized by the presence of O-linked N-acetylglucosamine moieties and a distinctive pentapeptide repeat (XFXFG). Nup153 is a peripheral NPC component that is implicated in protein and RNP transport and in the interaction of NPCs with the nuclear lamina. Nup153 contains a unique N-terminal region, a central domain consisting of four to five zinc fingers and a C-terminal region containing about 30 irregularly spaced FXFG repeats. Nup153 is cleaved by caspases during apoptosis. Nup153 interacts with TAP, which is essential for mRNA export, and associates with chromatin towards the end of anaphase, in parallel with the inner nuclear membrane protein LAP2. Nup153 is involved in NPC assembly, in anchoring NPCs within the nuclear envelope and in mediating specific nuclear import events.

# **Product:**

Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2

Mol	lecu	lar V	Veig	ht:	
154	kDa				

**Swiss-Prot:** 

P49790(Human)

**Purification&Purity:** 

ProA affinity purified

**Applications:** 

WB:1:500-1:2,000 ICC:1:100-1:500 FC:1:50-1:100


#### **Storage&Stability:**

Store at +4  $^{\circ}$ C after thawing. Aliquot store at -20  $^{\circ}$ C or -80  $^{\circ}$ C. Avoid repeated freeze / thaw cycles.


#### **Specificity:**

NUP153 polyclonal antibody detects endogenous levels of NUP153 protein.

# **DATA:**



Western blot analysis of NUP153 on K562 (1) and HepG2 (2) cell lysate using anti-NUP153 antibody at 1/500 dilution.



ICC staining NUP153 in A431 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

#### Note:

For research use only, not for use in diagnostic procedure.

# Bioworld Technology, Inc.

 Add:
 1660 South Highway 100, Suite 500 St. Louis Park, MN 55416,USA.

 Email:
 info@bioworlde.com

 Tel:
 6123263284

 Fax:
 6122933841

# Bioworld technology, co. Ltd.

 
 Add:
 No 9, weidi road Qixia District Nanjing, 210046, P. R. China.

 Email:
 info@biogot.com

 Tel:
 0086-025-68037686

 Fax:
 0086-025-68035151